


VARILITE Wheelchair Cushions

What is VARILITE ?


VARILITE cushions are an advanced range of self inflating air-foam flotation wheelchair cushions designed to offer comfort, support, durability and ease of use.

Overview of Technology

.. **Air-Foam Flotation**- refers to the combination of foam and air which ensures VARILITE can offer comfort, support and durability, all in a lightweight package, with most cushions weighing less than 1kg.



.. **Multi-stiffness foams**– most VARILITE cushions incorporate 3 stiffness's of foams; soft foam for under the coccyx and PSIS, medium foam for under the thighs and buttocks, and firm foam creating a supportive channel for the thighs.



.. **Self-Inflating**- All VARILITE cushions use self-inflating valves, making them easy to use and maintain. They require no pumps, no kneading, no manipulation, and cause no hassle.



.. **Cover**- The cushion cover provides extra comfort by encouraging heat, air, and vapour circulation. They are available in mesh or incontinence to ensure greater peace of mind.



So what

This means that VARILITE cushions do not have a downside; they don't compromise on comfort, support, durability, weight, maintenance, or pressure distribution. This not only makes them ideal for a wide range of wheelchair users, but also for every aspect of that wheelchair users diverse life. For instance, when going to the shops, their lightness makes them easy to transfer in and out of a car. When playing sport, the vibration dampening effect of air-foam floatation removes the "sea-sick sensation" experienced with many cushions. When simply sat the comfort and heat and vapour dissipation properties removes any risks or irritations. And regardless of what you get up to on your VARILITE, you can always be safe in the knowledge that you are not going to bottom out or needlessly risk skin ulcer development.

Without Compromise Range Highlights

Comfort
Pressure Distribution
Postural Support
Vibration Dampening
Low Maintenance
Temperature Transfer
Cushion Weight
Transfers
Cushion Choice
Cushion Life
Cushion Size
Incontinence Care
Ease of Use
User Weight
Warranty Period

VARILITE cushions come in a comprehensive range to meet all requirements.

.. **Pro-Form NX** offers a fully modifiable off the shelf seating solution, equipping seating specialists with an affordable tool to treat patients with limited hip flexion, obliquities, leg length discrepancies, or a tendency to scooter.



.. **Evolution** offers a high level of skin protection and support, which can be added to with the addition of a positioning base or wedge. The built in and modifiable wedge also ensures pressure is correctly distributed down the thighs, easing pressure on the delicate PSISs



.. **Zoid** is a comprehensive cushion for the sports user, building on the lightweight pedigree of the range. Designed to accommodate modified sport's required sitting position and ensure uninhibited movement, this is cleared to Olympic standards in most sports.



.. **Reflex** takes the plug and play idea behind *VARILITE* one step further, making it the ultimate throw and sit cushion. Providing advanced comfort to low risk users, the Reflex automatically adjusts to the correct air level for the each and every time.



VARILITE Wheelchair Cushions without Compromise



Health Care Innovations
AUSTRALIA

3/ 10-12 Carsten Rd, Gepps Cross, SA 5094.
Phone: 1300 499 282
Fax: 08 8125 5990
Email: sales@hcibiz.com.au
Web: www.healthcareinnovations.com.au