

Jay Cushions



Jay Easy

The Jay Easy is a lightweight cushion, designed for the active client who requires mild to moderate positioning and moderate to high skin protection. The pre-contoured foam provides a stable base of support for the pelvis. Medial and lateral thigh contouring promotes optimal leg positioning.



Jay Lite

The Jay Lite is the lightest cushion in its class constructed with 3DX™ Microclimatic spacer fabric and air flow base. It provides heat and moisture dissipation, good pressure distribution and stability. Designed for the client with low to moderate positioning and skin protection needs.



Jay Duo


The Duo lowers the wheelchair's seat to floor height and promotes foot propelling and / or independent transfers. Fluid in the rear protects ischials during foot propulsion.

Jay Cushions




Jay GS

The Jay GS is a growable paediatric cushion. It is designed for the child with high positioning and skin protection needs and features a structured positioning base with pressure relieving fluid overlay.



Jay Soft Combi P

Pre contoured foam base features a mild contouring for clients with low positioning / pressure care needs but high comfort needs. Incontinent cover comes standard.



Jay Basic

The Basic features durable molded foam providing mild contouring. The pre contoured foam base provides comfort and mild support. Incontinent cover comes standard.



Sunrise Medical Pty Ltd
Unit 10, 15 Carrington Rd
Castle Hill NSW 2154
Ph: 02 8858 9191
Fax - Admin: 02 8850 4711
Fax - Orders: 02 8858 9155
enquiries@sunrisemedical.com.au
www.sunrisemedical.com.au



JAY Backs & Cushions



Superior Clinical Seating.



Jay Backs



Jay Basic

Designed for a client with mild postural needs, this back provides mild posterior support and comfort. It is available in three different height options and either as a back replacement or to be used in front of wheelchair upholstery.



Jay J3

Offers a comprehensive back range to fully meet functional needs of wheelchair users featuring many sizes and easy attachment hardware. The Jay J3 is designed for moderate to high postural needs. Four support shapes offer full selection of contours to meet clients needs.



Jay J2

The Jay J2 is designed for moderate postural needs and facilitates thoracic extension. It has an angle adjustability of up to 15° and an optional spinal fluid pad insert that provides pressure and shear reduction.



Jay J2 Deep Contour

Designed for back support that provides both pelvic and trunk positioning. This back is for a client with moderate trunk and pelvic positioning needs.



Jay Focus Point

Designed for clients with high support needs, especially laterally. This back has a three piece contoured design that features a contoured aluminium shell and two multi-adjustable lateral wings. Scapular cut-outs minimise interference with upper extremity range of motion. It features tough durable 4 point mounting hardware.



Encompass Back


The custom version of a Focus Point Back. Customise with foam choices and lateral sizes to meet your client's exact specifications. Choose swing-away laterals for easier transfers.

Jay Cushions




Jay J3 with Fluid Pad

The J3 Cushion with fluid pad features outstanding pressure distribution, comfort, heat and moisture dissipation and optimal stability. It is comfortable, easy-to-order and can accommodate changing needs with customisable well shape and Jay fluid amounts. Jay J3 Fluid pad is appropriate for moderate to high risk clients that require stability and positioning in an easy to use low maintenance cushion.



Jay J3 with Air Pad

The J3 Cushion with air pad features outstanding pressure distribution, comfort, heat and moisture dissipation and optimal stability. It is comfortable, easy-to-order and can accommodate changing needs. Jay J3 Air pad is appropriate for moderate to high risk clients that require air insert for optimal pressure reduction but also benefit from stability offered by the Jay cushion base.



Jay J2 Deep Contour

Low maintenance cushion designed for the client at high risk for skin breakdown. The J2 deep contour has a deep contoured foam base which creates a foundation for postural stability and high level of fluid for maximal skin protection.



Jay J2

A pre contoured foam base creates a foundation for postural stability. The foam is closed cell and extremely lightweight. The base and fluid pad can be customised to correct or accommodate postures. Designed for clients with moderate to high risk.